


ACUPRIME®

EXCELLENCE AT SENSIBLE PRICES™



PRODUCT CATALOGUE

2025



About AcuPrime

AcuPrime is based in the historic city of Exeter, in the heart of the South West.

Since 2001, we have been aspiring to grow our company by providing our services to the professional community.

With the meridiuS® range, we strive to improve interconnectivity between the world and the body. We want to encourage people to live well.

We are the exclusive European distributor of DongBang® as well as the exclusive UK distributor of SpoBand® and Dr. Life intermittent pneumatic compression devices.

Our aim is to bring "Excellence at Sensible Prices" to healthcare professionals, helping restore vitality and wellbeing in both individuals and society.

www.acuprime.com

+44 (0)1392 829500

info@acuprime.com

01

Acupuncture Needles

21

Physiotherapy

25

Cupping Therapy & Massage

38

Press Pellet & Probes

48

Electrical Devices

53

Moxa & Heating Devices

69

Treatment Room Equipment

DongBang Acupuncture Needles

DongBang DB100



Silicone Coated Needles with Tubes & Spring Handles

- One needle with a guide tube per blister
- The needle features a spring handle and silicone coating
- Packed in 100 pcs per box
- EO gas sterilized - 5 year shelf life

Code	Size
DB1001615	0.16 x 15 mm
DB1001630	0.16 x 30 mm
DB1001815	0.18 x 15 mm
DB1001830	0.18 x 30 mm
DB1002015	0.20 x 15 mm
DB1002025	0.20 x 25 mm
DB1002030	0.20 x 30 mm
DB1002040	0.20 x 40 mm
DB1002050	0.20 x 50 mm
DB1002215	0.22 x 15 mm
DB1002230	0.22 x 30 mm
DB1002240	0.22 x 40 mm
DB1002525	0.25 x 25 mm
DB1002530	0.25 x 30 mm
DB1002540	0.25 x 40 mm
DB1002550	0.25 x 50 mm
DB1003030	0.30 x 30 mm
DB1003040	0.30 x 40 mm
DB1003050	0.30 x 50 mm
DB1003540	0.35 x 40 mm



DongBang DB102



Silicone Coated Needles without Guide Tubes

- One needle without a guide tube per blister
- The needle features a spring handle and silicone coating
- Packed in 100 pcs per box
- EO gas sterilized - 5 year shelf life

Code	Size
DB1021615	0.16 x 15 mm
DB1021630	0.16 x 30 mm
DB1021815	0.18 x 15 mm
DB1021830	0.18 x 30 mm
DB1022015	0.20 x 15 mm
DB1022025	0.20 x 25 mm
DB1022030	0.20 x 30 mm
DB1022040	0.20 x 40 mm
DB1022050	0.20 x 50 mm
DB1022230	0.22 x 30 mm
DB1022240	0.22 x 40 mm
DB1022515	0.25 x 15 mm
DB1022525	0.25 x 25 mm
DB1022530	0.25 x 30 mm
DB1022540	0.25 x 40 mm
DB1022550	0.25 x 50 mm
DB1023015	0.30 x 15 mm
DB1023030	0.30 x 30 mm
DB1023040	0.30 x 40 mm
DB1023050	0.30 x 50 mm
DB1023515	0.35 x 15 mm

DongBang Acupuncture Needles

DongBang DB105



Made from surgical stainless steel and EO gas sterilised, these sterile needles have a spring handle and easy-glide insertion action. The needles are gently released from the guide tube by removing the plastic tab.

- 5 needles and one guide tube per blister
- Disposable for single use
- 500pcs per box
- Silicon coated
- EO gas sterilized - 5 year shelf life

Code	Size
DB1051615	0.16 x 15 mm
DB1051630	0.16 x 30 mm
DB1051815	0.18 x 15 mm
DB1051830	0.18 x 30 mm
DB1052015	0.20 x 15 mm
DB1052025	0.20 x 25 mm
DB1052030	0.20 x 30 mm
DB1052040	0.20 x 40 mm
DB1052050	0.20 x 50 mm
DB1052215	0.22 x 15 mm
DB1052225	0.22 x 25 mm
DB1052230	0.22 x 30 mm
DB1052240	0.22 x 40 mm
DB1052515	0.25 x 15 mm
DB1052525	0.25 x 25 mm
DB1052530	0.25 x 30 mm
DB1052540	0.25 x 40 mm
DB1052550	0.25 x 50 mm
DB1053030	0.30 x 30 mm
DB1053040	0.30 x 40 mm
DB1053050	0.30 x 50 mm
DB1053060	0.30 x 60 mm



Acupuncture Needles

DongBang DB106



Made from surgical stainless steel and EO gas sterilised, these sterile needles have a spring handle and easy-glide insertion action. The needles are gently released from the guide tube by removing the plastic tab.

- 10 needles and one guide tube per blister
- Disposable for single use
- 1000pcs per box
- Silicon coated
- EO gas sterilized - 5 year shelf life



Code	Size
DB1061615	0.16 x 15 mm
DB1061630	0.16 x 30 mm
DB1061815	0.18 x 15 mm
DB1061830	0.18 x 30 mm
DB1061840	0.18 x 40 mm
DB1062015	0.20 x 15 mm
DB1062025	0.20 x 25 mm
DB1062030	0.20 x 30 mm
DB1062040	0.20 x 40 mm
DB1062050	0.20 x 50 mm
DB1052215	0.22 x 15 mm
DB1052225	0.22 x 25 mm
DB1062230	0.22 x 30 mm
DB1062240	0.22 x 40 mm
DB1062515	0.25 x 15 mm
DB1062525	0.25 x 25 mm
DB1062530	0.25 x 30 mm
DB1062540	0.25 x 40 mm
DB1062550	0.25 x 50 mm
DB1062560	0.25 x 60 mm
DB1063015	0.30 x 15 mm
DB1063030	0.30 x 30 mm
DB1063040	0.30 x 40 mm
DB1063050	0.30 x 50 mm
DB1063060	0.30 x 60 mm
DB1063530	0.35 x 30 mm
DB1063540	0.35 x 40 mm

DongBang Long Needles

- Long needle made from surgical stainless steel
- Ten spring handled needles with one tube per pouch



Code	Size	Contents
DB11175	0.35 x 75 mm	1000 pcs / box
DB11190	0.35 x 90 mm	1000 pcs / box
DB111105	0.35 x 105 mm	1000 pcs / box
DB111120	0.35 x 120 mm	1000 pcs / box
DB111135	0.35 x 135 mm	1000 pcs / box

DongBang Press Needles

- Press needles for ear and body acupuncture
- Made from surgical stainless steel



Code	Size	Contents
DB130S	0.20 x 1.0 mm	100 pcs / box
DB130M	0.20 x 1.5 mm	100 pcs / box

DongBang Intradermal Needles

- Silicone-free, stainless steel intradermal needles
- Packed in 50 pcs per box and 5 pcs per pouch



Code	Size
DB131S	0.12 x 3 mm
DB131M	0.12 x 5 mm
DB131L	0.14 x 7 mm

DongBang Hand Needles

- Hand acupuncture needles featuring spring handles
- Packed in 100 pcs per box and 10 pcs per pouch



Code	Size
DB132	0.18 x 8 mm

Acupuncture Needles

Meridius Acupuncture Needles

Meridius MCT Series Copper Handle Needles



- Chinese-style needle with copper wire handle
- One needle with a tube per blister
- Packed in 100 pcs per box

Code	Size
MCT2015	0.20 x 15 mm
MCT2213	0.22 x 13 mm
MCT2525	0.25 x 25 mm
MCT2530	0.25 x 30 mm
MCT2540	0.25 x 40 mm
MCT3030	0.30 x 30 mm
MCT3050	0.30 x 50 mm
MCT3075	0.30 x 75 mm
MCT3550	0.35 x 50 mm
MCT3575	0.35 x 75 mm

Meridius MST Series Silver Handle Needles



- Chinese-style needle with silver-plated wire handle
- One needle with a tube per blister
- Packed in 100 pcs per box

Code	Size
MST1615	0.16 x 15 mm
MST1630	0.16 x 30 mm
MST2015	0.20 x 15 mm
MST2525	0.25 x 25 mm
MST2530	0.25 x 30 mm
MST2540	0.25 x 40 mm
MST3040	0.30 x 40 mm
MST3060	0.30 x 60 mm
MST3075	0.30 x 75 mm
MST3575	0.35 x 75 mm

Meridius MC Series Copper Handle Needle

- Chinese-style needle with copper wire handle
- One needle without tube per blister
- Packed in 100 pcs per box



Code	Size
MC2015	0.20 x 15 mm
MC2525	0.25 x 25 mm
MC2540	0.25 x 40 mm
MC3030	0.30 x 30 mm

Meridius MS Series Silver Handle Needle

- One needle without guide tube per blister
- Packed in 100 pcs per box
- Features silicone coating for easy insertion



Code	Size
MS2015	0.20 x 15 mm
MS2525	0.25 x 25 mm
MS2530	0.25 x 30 mm
MS2540	0.25 x 40 mm
MS3030	0.30 x 30 mm
MS3040	0.30 x 40 mm

Meridius ME130 Press Needles

- Press needles for ear and body acupuncture
- Made from surgical stainless steel



Code	Size	Contents
ME130S2010	0.20 x 1.0 mm	100 pcs / box
ME130M2015	0.20 x 1.5 mm	100 pcs / box

ASP Original Auriculotherapy Needles



The ASP Original needle is composed of a small stainless steel needle housed in a plastic injector. This injector allows the practitioner to easily and precisely implant the needle into the patient's ear. Inserting the needle in this way avoids any risk of contamination. After placement, the needle can be protected by an adhesive plaster. The needle ensures a permanent action for as long as it remains in place. A black magnet housed in the injector can be turned around in front of the implanted needle. ASP Original needles are gamma ray sterilized and for single use only.

Invented in collaboration with Dr Paul Nogier, the ASP Original or "Aiguille Semi-Permanente" is used to perform ear acupuncture. It has been marketed for more than 30 years by Sedatelec, and the shape of this small needle makes it exceptionally efficient.

Code	Size
ASC8	Classic 8 pcs/box
ASC80	Classic 80 pcs/box
ASC200	Classic 200 pcs/box
ASG8	Gold 8 pcs/box
ASG80	Gold 80 pcs/box
ASG200	Gold 200 pcs/box

ORCA™ Acupuncture Needles

ORCA™ Accent Copper Handled Needles

- Chinese-style needle with copper wire handle
- One needle with a tube per blister
- Packed in 100 pcs per box



Code	Size
OCT1630	0.16x30 mm
OCT2015	0.20x15 mm
OCT2030	0.20x30 mm
OCT2213	0.22x13 mm
OCT2225	0.22x25 mm
OCT2525	0.25x25 mm
OCT2530	0.25x30 mm
OCT2540	0.25x40 mm
OCT3013	0.30x13 mm
OCT3025	0.30x25 mm
OCT3030	0.30x30 mm
OCT3040	0.30x40 mm
OCT3050	0.30x50 mm
OCT3075	0.30x75 mm
OCT3090	0.30x90 mm
OCT3525	0.35x25 mm
OCT3530	0.35x30 mm
OCT3540	0.35x40 mm
OCT3550	0.35x50 mm
OCT3560	0.35x60 mm
OCT3570	0.35x70 mm
OCT3575	0.35x75 mm
OCT35100	0.35x100 mm

ORCA™ Advanced Copper Handled Needles

- Chinese-style needle with copper wire handle
- One needle without tube per blister
- Packed in 100 pcs per box



Code	Size
OC2015	0.20x15 mm
OC2515	0.25x15 mm
OC2525	0.25x25 mm
OC2530	0.25x30 mm
OC2540	0.25x40 mm
OC25100	0.25x100 mm
OC3030	0.30x30 mm
OC3050	0.30x50 mm

ORCA™ Accent Silver Handled Needles

- Chinese-style needle with silver-plated wire handle
- One needle with a tube per blister
- Packed in 100 pcs per box



Code	Size
OST1630	0.16x30 mm
OST1830	0.18x30 mm
OST1840	0.18x40 mm
OST2015	0.20x15 mm
OST2213	0.22x13 mm
OST2515	0.25x15 mm
OST2525	0.25x25 mm
OST2530	0.25x30 mm
OST2540	0.25x40 mm
OST2570	0.25x70 mm
OST3030	0.30x30 mm
OST3040	0.30x40 mm
OST3050	0.30x50 mm
OST3060	0.30x60 mm
OST3070	0.30x70 mm
OST3075	0.30x75 mm
OST3080	0.30x80 mm
OST3090	0.30x90 mm
OST30100	0.30x100 mm
OST30120	0.30x120 mm
OST3275	0.32x75 mm
OST3530	0.35x30 mm
OST3540	0.35x40 mm
OST3550	0.35x50 mm
OST3560	0.35x60 mm
OST3575	0.35x75 mm
OST35100	0.35x100 mm
OST35120	0.35x120 mm
OST35130	0.35x130 mm
OST35135	0.35x135 mm

ORCA™ Advanced Silver Handled Needles



Chinese-style needles with silver-plated wire handles

- One needle without guide tube per blister
- Packed in 100 pcs per box
- Features silicone coating for easy insertion
- EO gas sterilized

Code	Size
OS1615	0.16x15 mm
OS2015	0.20x15 mm
OS2515	0.25x15 mm
OS2525	0.25x25 mm
OS2530	0.25x30 mm
OS2540	0.25x40 mm
OS3030	0.30x30 mm
OS3040	0.30x40 mm
OS3050	0.30x50 mm

ORCA™ Premium Plastic Handled Needles



Coloured-coded plastic handled needles

- One needle without tube per blister
- Packed in 100 pcs per box

Code	Size
OP1215	0.12x15 mm
OP1415	0.14x15 mm
OP1607	0.16x07 mm
OP1613	0.16x13 mm
OP1630	0.16x30 mm
OP1815	0.18x15 mm
OP2007	0.20x07 mm
OP2013	0.20x13 mm
OP2015	0.20x15 mm
OP2030	0.20x30 mm
OP2207	0.22x07 mm
OP2213	0.22x13 mm
OP2525	0.25x25 mm
OP2540	0.25x40 mm
OP3030	0.30x30 mm
OP3050	0.30x50 mm

Acupuncture Needles

ORCA™ Acupuncture Needles

ORCA™ Resolve Silver Handled Needles



Chinese-style needles with silver-plated wire handles

- Five needles with one guide tube per blister
- Packed in 500 pcs per box
- Features silicone coating for easy insertion
- EO gas sterilized

Code	Size
OS52525	0.25x25 mm
OS52530	0.25x30 mm
OS52540	0.25x40 mm
OS52550	0.25x50 mm
OS53030	0.30x30 mm
OS53040	0.30x40 mm

ORCA™ Resolve Copper Handled Needles



Chinese-style needles with copper wire handles

- Five needles with one guide tube per blister
- Packed in 500 pcs per box
- Features silicone coating for easy insertion
- EO gas sterilized

Code	Size
OC52525	0.25 x 25 mm
OC52530	0.25 x 30 mm

ORCA™ Natural Copper Handled Needles



Silicone-free needles with copper wire handles

- One needle without guide tube per blister
- Packed in 100 pcs per box
- Ideal for electro-acupuncture
- EO gas sterilized

Code	Size
ON2015	0.20x15 mm
ON2213	0.22x13 mm
ON2225	0.22x25 mm
ON2525	0.25x25 mm
ON2530	0.25x30 mm
ON2540	0.25x40 mm
ON3030	0.30x30 mm
ON3050	0.30x50 mm

Seirin® Acupuncture Needles

Seirin® J-Type Needles



- One needle with a guide tube per blister
- Packed in 100 pcs per box
- Every needle features a skin-friendly plastic handle
- EO gas sterilized

Code	Size
SJ1215	0.12x15 mm
SJ1230	0.12x30 mm
SJ1415	0.14x15 mm
SJ1430	0.14x30 mm
SJ1615	0.16x15 mm
SJ1630	0.16x30 mm
SJ1640	0.16x40 mm
SJ1830	0.18x30 mm
SJ1840	0.18x40 mm
SJ2030	0.20x30 mm
SJ2050	0.20x50 mm
SJ2530	0.25x30 mm
SJ2540	0.25x40 mm
SJ2550	0.25x50 mm
SJ3030	0.30x30 mm
SJ3040	0.30x40 mm
SJ3050	0.30x50 mm
SJ3060	0.30x60 mm



Seirin® Spinex Intradermal Needles



- One EO gas sterilized needle per blister
- Packed in 100 pcs per box
- Features a ring handle measuring 1.5 mm in diameter
- Made from surgical stainless steel

Code	Size
DBSSP1203	0.12 x 3 mm
DBSSP1204	0.12 x 4 mm
DBSSP1406	0.14 x 6 mm

Seirin® B-Type Needles



- One EO gas sterilized needle per blister
- Packed in 100 pcs per box
- Features silicone coating for easy insertion
- Made from surgical stainless steel
- Every B-Type needle features a lightweight plastic handle

Code	Size
SB1615	0.16 x 15 mm
SB2015	0.20 x 15 mm
SB2540	0.25 x 40 mm
SB3030	0.30 x 30 mm

Seirin® New Pyonex Needles



The patented plastic packaging allow the practitioner to insert the needle without touching the needle tip or the surface of the Micropore™ adhesive tape.

The skin friendly adhesive allows the skin to breathe during treatment. With Seirin's smooth insertion, the result is a painless experience. Packed in 100 pcs / box.

Code	Size	Colour
DBSNP2003	0.20 x 0.3 mm	Orange
DBSNP2006	0.20 x 0.6 mm	Yellow
DBSNP2009	0.20 x 0.9 mm	Green
DBSNP2012	0.20 x 1.2 mm	Blue
DBSNP2015	0.20 x 1.5 mm	Pink

Acupuncture Accessories

Hand Needle Injectors



- Assists in the insertion of hand needles
- Magnets inside pick up, hold and secure the needle
- The weight of the plunger will drop, using gravity or the spring force to insert the needle
- DB135 allows adjustments for the depth of insertion

Code	Type
DB135	Spring Force

NIDO 7-Star Needles



- Disposable 7-star needle
- Pipe type double-ended – ideal for bleeding acupuncture
- Individually packed 7-star dermal needle
- Double-side needle-18cm(7") long with 7 needles

Code	Type
ND151HQ	1 Head

Acupuncture Needles

Lancet Devices

DongBang NANOLET™ Lancets



Code	Size
DB40023G	0.60 mm (23G)
DB40026G	0.45 mm (26G)
DB40028G	0.35 mm (28G)
DB40030G	0.30 mm (30G)

LANZO



Adjustable tip - 5 settings for custom lancet penetration
Lancet ejector - no need to touch the lancet to remove it
Reduced pain - thanks to the linear tracking system
Silent - the carrier is provided with two springs and a stopper

Optional clear tip for alternative site testing designed in order to reduce pain and secure a blood sample.

Code
DB403A

DongBang Lancet Device



- To be used with DB400 NANOLET™ Lancets
- Press the trigger on the top of the pen to push the needle

Code
DB403B

Park Sham Device

Park Sham Research Tools



Fig 1


The Park Sham Device is an apparatus used to administer the sham needle and real needles. The Park Sham Needle appears similar to the real needle, but is blunt and free to slide within the handle so that when pressed, it telescopes into the handle instead of penetrating the skin. The device holds the un-inserted sham needle in place and maintains the sterility of a genuine acupuncture needle.

This device consists of a guide tube: a special oversized tube with base flange which is equipped with double-sided tape and can therefore stick on target skin (Fig 1). Since the Park Sham Device was first developed in 1999, several teams of researchers in UK, Japan, Korea and USA have tested its credibility and validity in various settings.

Based on these findings over the past 19 years it has been frequently used in research projects, often funded by government agencies including the US National Institute of Health. Park and his fellow colleagues, who developed this device, state that this device is indistinguishable from the same procedure using real needles in acupuncture naive subjects and inactive, where the specific needle sensation is taken as a surrogate measure of activity. However, due to the nature of minimal physiological activity required to achieve successful subject blinding in acupuncture trial, they emphatically call for caution in application of sham acupuncture control and in interpretation of the findings. For this reason, Park is keen to call this a bona fide version of sham control rather than final version.

- Park J. et al. Development of a new sham needle, *Acupunct Med* 17 (1999), pp. 110-112.
- Park J. et al., Validating a new nonpenetrating sham acupuncture device: two randomised controlled trials, *Acupunct Med* 20 (2002), pp. 168-174.
- Whale CA et al. Pilot study to assess the credibility of acupuncture in acute exacerbations of chronic obstructive pulmonary disease. *Acupunct Med.* 2009 Mar;27(1):13-5.
- Tsukayama H et al. Factors that influence the applicability of sham needle in acupuncture trials: two randomized, single-blind, crossover trials with acupuncture-experienced subjects. *Clin J Pain.* 2006 May;22(4):346-9.
- Chae Yet al. The neural substrates of verum acupuncture compared to non-penetrating placebo needle: an fMRI study. *Neurosci Lett.* 2009 Jan 30;450(2):80-4.
- Tan CW et al. Discrimination of real and sham acupuncture needles using the Park sham device: a preliminary study, *Arch Phys Med Rehabil* 90 (2009), pp. 2141-2145.
- Downs NM et al. The effect of real and sham acupuncture on thermal sensation and thermal pain thresholds. *Arch Phys Med Rehabil* 2005;86:1252-7.
- Park J, et al. Acupuncture for subacute stroke rehabilitation: a Sham-controlled, subject and assessor blind, randomized trial. *Arch Intern Med.* 2005 Sep 26;165(17):2026-31.
- Smith Pet al. The efficacy of acupuncture in the treatment of temporomandibular joint myofascial pain: a randomised controlled trial. *J Dent.* 2007 Mar;35(3):259-67.
- Kennedy Set al. Acupuncture for acute non-specific low back pain: a pilot randomised non-penetrating sham controlled trial. *Complement Ther Med.* 2008 Jun;16(3):139-46.
- Park J. Sham needle control needs careful approach. *Pain.* 2004 May;109(1-2): 195-6
- Park JJ. Developing and validating a sham acupuncture needle. *Acupunct Med.* 2009 Sep;27 (3):93.

How to use the Park Sham Device



To avoid compromising the sterility, you could peel off only the top of the package to expose the handle of the needle.

Peel off the cover of the double-sided sticky tape and place the device on the skin while keeping the needle handle projecting about 5mm from the top of the Park tube.

Then tap the needle gently to insert it about 5mm followed by moving down the upper telescopic guide tube to allow sufficient exposure of the handle for manipulation of the needle.

To remove: first take out the needle and dispose of it safely, then peel the flange from the skin, with the guide tube in place. Doing this quickly will minimise discomfort.

Code	Product
PSD01	Park Sham Device
PSD03	Sham Needle Only
PSD04	Sticky Tape
PSD05	Tube Only
PSD06	Real Matching Needle

Needles Sizes summary chart

DongBang DB100

Length (mm)	Width (mm)						
	0.16	0.18	0.20	0.22	0.25	0.30	0.35
15							
25							
30							
40							
50							

DongBang DB102

Length (mm)	Width (mm)						
	0.16	0.18	0.20	0.22	0.25	0.30	0.35
15							
25							
30							
40							
50							
60							

DongBang DB105

Length (mm)	Width (mm)						
	0.16	0.18	0.20	0.22	0.25	0.30	0.35
15							
25							
30							
40							
50							
60							

DongBang DB106

Length (mm)	Width (mm)						
	0.16	0.18	0.20	0.22	0.25	0.30	0.35
15							
25							
30							
40							
50							
60							

Meridius MCT

Length (mm)	Width (mm)				
	0.20	0.22	0.25	0.30	0.35
13					
15					
25					
30					
40					
50					
75					

Meridius MST

Length (mm)	Width (mm)				
	0.16	0.20	0.25	0.30	0.35
15					
25					
30					
40					
60					
75					

Meridius MC

Length (mm)	Width (mm)		
	0.20	0.25	0.30
15			
25			
30			
40			

Meridius MS

Length (mm)	Width (mm)		
	0.20	0.25	0.30
15			
25			
30			
40			

Acupuncture Needles

Needles Sizes summary chart

ORCA Accent Copper

Length (mm)	Width (mm)					
	0.16	0.20	0.22	0.25	0.30	0.35
13						
15						
25						
30						
40						
50						
60						
75						
90						
100						

ORCA Accent Silver

Length (mm)	Width (mm)							
	0.16	0.18	0.20	0.22	0.25	0.30	0.32	0.35
13								
15								
25								
30								
40								
50								
60								
70								
75								
80								
90								
100								
120								
130								
135								

ORCA Premium

Length (mm)	Width (mm)							
	0.12	0.14	0.16	0.18	0.20	0.22	0.25	0.30
7								
13								
15								
25								
30								
40								
50								

Seirin J Type

Length (mm)	Width (mm)						
	0.12	0.14	0.16	0.18	0.20	0.25	0.30
15							
25							
30							
40							
50							
60							

Seirin B Type

Length (mm)	Width (mm)			
	0.16	0.20	0.25	0.30
15				
30				
40				

Akumata Standard – acupressure mat with pillow



- This mat helps alleviate muscle tension, improve circulation, and provide relief for back, neck, and joint discomfort.
- Made from 100% cotton and durable HIPS material, ensuring comfort and longevity.
- Length 73cm x Width 43cm x Thickness 2cm

Code	Colour
AM03BK-GD	Gold
AM03BK-GR	Green
AM03BK-MT	Mint

Akumata Long – acupressure mat with pillow



- This mat helps alleviate muscle tension, improve circulation, and provide relief for back, neck, and joint discomfort.
- Longer surface area, ensuring a more immersive experience for full-body relaxation
- Length 110cm x Width 43cm x Thickness 2cm

Code	Colour
AM02BK-GD	Gold
AM02BK-GR	Grey
AM02BK-MT	Mint

Akumata ECO – Ecological acupressure mat with pillow



- Relieves muscle and joint pain
- Natural linen and cotton, hypoallergenic HIPS materials
- Length 73cm x Width 43cm x Thickness 2cm

Code	Colour
AM01LN-BK	Black
AM01LN-GD	Gold

Kinesiology Tape Working Principles



Kinesiology taping was developed in the 1970s by Japanese chiropractor Kenzo Kase, who went to the U.S. in the mid-1990s. Unlike standard athletic taping, which often involves wrapping a joint for support and compression, kinesiology tape is placed in a variety of patterns depending on the injury. It is pulled to differing degrees of tension to create the desired effect and is typically worn for two to five days, unlike standard tape, which is used mainly during an activity.

Kinesiology tape can be used to pull back a shoulder that is hunching forward or to reduce swelling in a joint. Tape can be used to pull up the skin and create an area of low pressure where fluid can move and drain. Some clinicians believe the tape provides stimulation to skin cells that affects pain pathways - similar to rubbing a spot that hurts. (Summarized from the Wall Street Journal Apr 27, 2010)

Directions for use



Always shape the cut ends of the tape in such a way that they are curved with no sharp corners. This helps to reduce the problem of the corners catching in clothing, causing the tape to become dislodged.

To correctly apply the tape, first measure the required length while the backing tape is still in place, then shape the corners as outlined above. Position the area to be taped so that the underlying muscles are fully stretched. Remove the backing and lay the tape on the skin sticking it in place. Do not stretch the tape, simply lay it on the skin.

The exact length of time it remains in place will depend on how greasy the skin is and how much the subject perspires. It is possible to shower carefully, keeping the tape in place. The area should be patted dry and not rubbed. If the tape then feels unpleasantly damp, the area may be dried using a hairdryer on the lowest heat setting, held some distance from the body.

Kinesiology Tapes

meridius® ADAPTA™ Tape



Meridius® Adapta Tape is the result of years of development and testing conducted by top athletes and renowned scientists. This kinesiology tape is comfortable to wear, but also sticks for up to 4 days due to the premium hypo-allergenic adhesive that's being used.

Using only high quality materials, we have created a water resistant, bi-linearly elastic tape that doesn't compromise the natural range of motion. Adapta Tape has shown results in improving circulation, helping with muscle soreness, supporting joints and improving posture.



Code	Size	Colour
ME01BE	5 m	Desert Beige
ME03BK	5 m	Chaos Black
ME05BL	5 m	Arctic Blue
ME07PK	5 m	Toxic Pink
ME09YL	5 m	Atomic Yellow
ME11GR	5 m	Viper Green
ME12OR	5 m	Tiger Orange
ME13RD	5 m	Volcano Red
ME14VT	5 m	Mystic Violet

Code	Size	Colour
ME32BE	32 m	Desert Beige
ME32BK	32 m	Chaos Black
ME32BL	32 m	Arctic Blue
ME32PK	32 m	Toxic Pink
ME32YL	32 m	Atomic Yellow

AcuTape™



This tape is made of:

- Cotton 72%
- Spandex 3%
- Ethyl acetate 15%

Code	Size	Colour
DAT01BE	5 m	Beige
DAT03BK	5 m	Black
DAT05BL	5 m	Blue
DAT07PK	5 m	Pink
DAT09YL	5 m	Yellow
DAT11GR	5 m	Green

Resistance Bands

SpoBand® Elastic Resistance Band



SpoBand® is a premium resistance band that is used for muscle toning, posture correction, cardio exercises, rehabilitation and yoga.

While the typical rubber or latex bands can begin to smell after usage, SpoBand® doesn't. It's also easy to clean as it is machine washable.

Another important feature is that SpoBand® doesn't stick and melt as the regular bands do. Moreover, the high quality materials will not tear or snap like latex bands are prone to.

SpoBand features 9 loops for different resistances making it versatile, safe and easy to use.

Code	Name	Size
MESB01	SpoBand 20	110 x 2.0 x 0.5 cm
MESB02	SpoBand 25	110 x 2.5 x 0.5 cm
MESB03	SpoBand 30	110 x 3.0 x 0.5 cm
MESB04	SpoBand 35	110 x 3.5 x 0.5 cm
MESB05	SpoBand 45	110 x 4.5 x 0.5 cm
MESB06	SpoBand 55	110 x 5.5 x 0.5 cm

Cupping Therapy & Massage

meridius® Cupping Sets

meridius® 17pcs Cupping Set



meridius® 10pcs Cupping Set



Meridius® Cupping Sets are presented in three sizes: 10 pieces, 17 pieces, and 30 pieces. They come blue case with a firm handle which is ideal for travelling and storage. The blue silicon valve can be detached from the cup for easy sterilisation and maintenance.

Our cupping sets are extremely versatile due to the large number of cupping shapes and dimensions available. Made from high-quality polycarbonate, Meridius® Cupping Sets are extremely durable and do not crack from falling or prolonged usage.

ME17 Contents

Code	Diameter	Pieces
ME17-R01	Round, 55 mm	8
ME17-C02	Contoured, 43 mm	1
ME17-C03	Contoured, 41 mm	1
ME17-C04	Contoured, 34 mm	1
ME17-R04	Round, 27 mm	1
ME17-TR	Triangle, 23 mm	1
ME17-OV1	Contoured, 72 mm	1
ME17-C01	Contoured, 61 mm	1
ME17-L01	Round, 50 mm	1
ME17-L02	Round, 40 mm	1
ME17-T	Extension Tube	1
ME17-P	Hand Pump	1

ME10 Contents

Code	Diameter	Pieces
ME17-R01	Round, 55 mm	6
ME17-C02	Contoured, 43 mm	2
ME17-C01	Contoured, 61 mm	1
ME17-C04	Contoured, 34 mm	1
ME17-P	Hand Pump	1
ME17-T	Extension Tube	1

meridius® 30pcs Cupping Set



ME30 Contents

Code	Diameter	Pieces
ME17-C01	Contoured, 61 mm	2
ME17-C02	Contoured, 43 mm	2
ME17-C03	Contoured, 41 mm	2
ME17-C04	Contoured, 34 mm	2
ME17-R04	Round, 27 mm	1
ME17-TR	Triangle, 23 mm	1
-	Square, 6 x 4 cm	20
ME17-T	Extension Tube	1
ME17-P	Hand Pump	1

Cupping Therapy & Massage

Replacement cups



ME17-R01



ME17-C02



ME17-C03



ME17-C04



ME17-R04



ME17-TR



ME17-OV1



ME17-C01



ME17-L01



ME17-L02



ME17T



ME17P

Code	Shape	Diameter
ME17-R01	Round	Round, 55 mm
ME17-C02	Contoured	Contoured, 43 mm
ME17-C03	Contoured	Contoured, 41 mm
ME17-C04	Contoured	Contoured, 34 mm
ME17-R04	Round	Round, 27 mm
ME17-TR	Triangle	Triangle, 23 mm
ME17-OV1	Contoured	Contoured, 72 mm
ME17-C01	Contoured	Contoured, 61 mm
ME17-L01	Round	Round, 50 mm
ME17-L02	Round	Round, 40 mm
ME17-T	Extension Tube	Standard
ME17-P	Hand Pump	Standard

meridius® VIVA Cupping Set, 10 & 16 cups



VIVA Cupping Set Contents

Code	Diameter
ME-VIVA10	Round, 58mm x 2 pcs
	Round, 48mm x 4 pcs
	Round, 36mm x 1 pcs
	Round, 26mm x 1 pcs
	Round, 20mm x 1 pcs
	Round, 16mm x 1 pcs
	Hand Pump
ME-VIVA16	Round, 58mm x 4 pcs
	Round, 48mm x 8 pcs
	Round, 36mm x 1 pcs
	Round, 26mm x 1 pcs
	Round, 20mm x 1 pcs
	Round, 16mm x 1 pcs
	Tube + Hand Pump

- Available in 10 pieces set and 16 pieces set.
- New designed detachable silicone valve for easy cleaning and sterilization.
- Durable and long lasting cups made of Polycarbonate.

Replacement cups



ME-VIVA-01



ME-VIVA-02



ME-VIVA-03



ME-VIVA-04



ME-VIVA-05



ME-VIVA-06

Code	Shape	Diameter
ME-VIVA-01	Round	58mm
ME-VIVA-02	Round	48mm
ME-VIVA-03	Round	36mm
ME-VIVA-04	Round	26mm
ME-VIVA-05	Round	20mm
ME-VIVA-06	Round	16mm

Cupping Therapy & Massage

DongBang Cupping Sets

DongBang Plastic Cupping Set



- Elegant lightweight case - Easy to carry and store
- Cups are non-breakable with valve on the top
- Place the pump on top of the cup and apply pressure
- To release pressure press the valve

Code	Diameter
DB350	51 mm x 13 pcs
	46 mm x 1 pcs
	41 mm x 1 pcs
	33 mm x 1 pcs
	25 mm x 1 pcs
	Tube + Hand Pump
DB360	51 mm x 6 pcs
	46 mm x 1 pcs
	41 mm x 1 pcs
	33 mm x 1 pcs
	25 mm x 1 pcs
	Hand Pump

DongBang Disposable Cups



- Disposable plastic cups for single use
- Compatible with DB304 hand pistol pump

Code	Diameter	Contents
DB310-1	55 mm	5 cups / pack
DB310-2	51 mm	5 cups / pack
DB310-3	46 mm	5 cups / pack
DB310-4	41 mm	5 cups / pack
DB310-5	36 mm	5 cups / pack

Bio Magnet Cupping Set



- A 17 piece plastic cupping set with magnets which also apply a mild acu-pressure to the body
- Each set includes 7 red and 8 blue colour coded, 0.3" dia, 1,000 gauss, spring adjustable, 0.3" dia x 1.7" long, removable magnetic insert probes. These fit the 2" and 1.5" dia cups
- Also includes pistol grip hand pump, extension tube and carrying case

Code
DBA341

Plastic Cupping Set



Code	Size
DBA340-6	65 mm x 1 pcs
	55 mm x 1 pcs
	45 mm x 1 pcs
	35 mm x 1 pcs
	Tube + Pump
DBA340-12	65 mm x 2 pcs
	55 mm x 2 pcs
	45 mm x 2 pcs
	35 mm x 4 pcs
	Tube + Pump
DBA340-24	65 mm x 12 pcs
	55 mm x 4 pcs
	45 mm x 4 pcs
	35 mm x 4 pcs
	Oil + Tube + Pump

Cupping Therapy & Massage

Glass Cupping Sets

DaeGun Glass Cupping Set



- Silver colour carrying case with prime quality glass cups
- Offers better cupping treatment
- Glass cup with valve on the top
- Contains 14 cups with a hand pump

Code
DB302

Glass Fire Cupping Jars



- Thick, heavy glass jar for cupping therapy
- Ideal for fire cupping
- Burn alcohol cotton balls to create a vacuum
- Easy to clean

Code	Diameter
DB305-XL	75 mm
DB305-L	65 mm
DB305-M	55 mm
DB305-S	45 mm
DB305-XS	40 mm

Other Cupping Sets

HACI® Magnetic Cupping Set



- Rubber bulb plastic suction cups
- 2300 gauss strength magnets
- 5 bio-north cups and 5 bio-south cups

Code
DBA320

Silicone Cupping Jars



- Each set contains 4 cups (25, 40, 50, 60mm diameter)

Code	Colour
DBA322T	Transparent
DBA322	Blue

Bamboo Cupping Jars



- Natural bamboo cups in three different sizes
- Cups can be soaked in water for infusion with herbs

Code	Contents
DBA321	3 pcs / box

Silicone Cups for Face and Eyes



Code	Size	Colour
DBA325S-BL	6 mm x 24 mm	Blue
DBA325S-GR	6 mm x 24 mm	Green
DBA325S-PK	6 mm x 24 mm	Pink
DBA325S-PL	6 mm x 24 mm	Purple
DBA325S-CL	6 mm x 24 mm	Clear

Silicone Cups for General Use



Code	Size	Colour
DBA324L-01BL	47 mm x 80 mm	Blue
DBA324L-01CL	47 mm x 80 mm	Clear
DBA324L-01GR	47 mm x 80 mm	Green
DBA324L-01PK	47 mm x 80 mm	Pink
DBA324L-01PL	47 mm x 80 mm	Purple
DBA324M-01BL	33 mm x 60 mm	Blue
DBA324M-01CL	33 mm x 60 mm	Clear
DBA324M-01GR	33 mm x 60 mm	Green
DBA324M-01PK	33 mm x 60 mm	Pink
DBA324M-01PL	33 mm x 60 mm	Purple
DBA324S-01BL	19 mm x 80 mm	Blue
DBA324S-01CL	19 mm x 80 mm	Clear
DBA324S-01GR	19 mm x 80 mm	Green
DBA324S-01PK	19 mm x 80 mm	Pink
DBA324S-01PL	19 mm x 80 mm	Purple

Medical grade silicone, odour free and hygienic. Transparency enables observation of the skin and soft tissue inside the cups during the application.

How do they work?

Expelling air out of the cup creates negative pressure and causes the skin and soft tissue to rise up. Thanks to the soft rims of the silicone cups, these cups can be applied to contoured and other tricky areas of the body. The soft rim creates a flawless seal around the connected skin and is far more ergonomic than its glass counterparts.

How are they used?

Ensure the skin is clean. Apply massage cream, gel or oil to moisten the skin and prepare it for both massaging and cup gliding. Squeeze the cup in the middle to release air from inside. Place the squeezed cup onto the skin and hold it against the skin to ensure a complete seal is created.

Slowly release and allow the pressure to lift the skin and flesh. Grasp the base of the cup and gently massage or glide along the skin to create a wave of lifting. Massage oil, cream, gel, or soapy water can be considered to keep the skin lubricated to assist the massaging.

Non-Sterile Cotton Wool Balls



Code	Size
DBA6059000	Small, 500 pcs / pack
DBA6059001	Large, 250 pcs / pack

Sterile Cotton Wool Balls



Code	Size
DBA6059008	5 pcs / pack

Rectangle W hite Jade Gua Sha Board



Code	Size
BN6011CW	5.5 x 8 cm

Comb Shaped Rectangle Gua Sha Tool



Code	Size
BN6030-05	5 pcs / pack

Gua Sha Board



- Flat rectangle with slight contour at sides for traditional Gua Sha treatments

Code	Size
BN6007	3.5 x 6.5 cm
BN6008	4.5 x 8 cm
BN6009	5.5 x 10 cm

Gua Sha Jade



Code	Size
BN6011H	5 x 7.5 cm
BN6011C	5.5 x 8 cm
BN6011F	4 x 11 cm

Quartz Gua Sha Rod



Code
BN6012

Arched Fish Shaped Gua Sha Board



Code	Size
BN6025-FISH	15 x 5 cm

Gua Sha Board Shell Shape



Code	Size
BN6016	8 x 5 cm

Horn shape 15 cm long with one flat



Code	Size
BN6030-03	15 cm x 1 cm ø

Gua Sha Roller Jade



Code
BN6015D

Gua Sha Agate Facial Roller



Code
BN6016

Gua Sha Massage Tool With Finger Hole



Code	Size
BN6030-01	5 x 3 cm

Curved Horn Shape Gua Sha Rod



Code	Size
BN6001-160	5 x 3 cm

Press Pellet

Magnetic Ball Plasters



- 800 gauss strength
- 2mm diameter with tan-coloured surgical tape backing

Code	Type	Contents
DBA424	Silver plated, 2mm	100pcs / box
DBA425	Gold plated, 2mm	100pcs / box

Magrain Ion Pellet



- Sakamura Magrain are used for both ear and body points
- Magrain are silver colour 1.2mm alloy granules
- Magrain is unique in that it provides an ionizing reaction which further enhances the compression stimulation to the point

Code	Type	Contents
DB509006C	Non-plated	300pcs / pack
DB509006GC	Gold plated	300pcs / pack

Vaccaria Ear Seeds



- Vaccaria Ear Seeds are attached to an individual plaster, so that after the removal they can be placed immediately
- Proper application on the desired place is possible thanks to the plaster
- Typical application of the ear seed is on the acupuncture points of the ear, in order to achieve the lasting effect of the pressure
- Based on their size they can also be attached to other locations of the body

Code	Size	Contents
DBA426	2 mm	100pcs / box

Press Pellet



- Ideal for Korean Hand Acupuncture Therapy
- The metal protrusions provide gentle stimulation
- Pellets may remain in place for several hours

Code	Type	Contents
DB506A	Silver plated	100pcs / box
DB506B	Silver plated	30pcs / box
DB506C	Gold plated	100pcs / box
DB506D	Gold plated	30pcs / box

Tourmaline Stone Granules for Ear



- The product must be used by a qualified practitioner
- Product is for single use only, with care, to attach securely to ear
- Store in a safe place away from children, chemicals, moisture and extreme temperature

Code	Type	Contents
DBA531	1 ball	420pcs / box
DBA531-123	1 ball	126pcs / box
	2 balls	84pcs / box
	3 balls	140pcs / box

Magnets

HANIL Magnetic Press Pellet



- Magnets for acupuncture application
- 5mm diameter pellet on plaster backing
- 800 gauss
- Contains 100 pcs per box

Code	Diameter	Contents
DB511	5 mm	100pcs / box

ACCU-BAND 800G



ACCU-BAND gives a constant pressure and a concentrated magnetic flux of 800 gauss, 6000 gauss, 9000 gauss. Adheres even in water and will not rust. The south pole faces towards the body.

Code	Type
DB509001G	Gold plated, 800 gauss, 24pcs/ box
DB509001S	Non-plated, 800 gauss, 24pcs/ box
DB509002S	Gold plated, 6000 gauss, 12pcs/ box
DB509003S	Gold plated, 9000 gauss, 12pcs/ box

Hpro Ferrite Magnets



Code	Type
DBA5901	Non-plate, 700 gauss 10pcs/Box
DBA5901G	Gold-plated, 700 gauss 10pcs/Box
DBA5902	Non-plate, 1000 gauss 10pcs/Box
DBA5902G	Gold-plated, 1000 gauss 10pcs/Box
DBA5905	Non-plated, 1200 gauss 10pcs/Box
DBA5905G	Gold-plated, 1200 gauss 10pcs/Box
DB5903G	Gold-plated, 2500 gauss 10pcs/Box

Meridius AcuMagnet Patch

- **Effective Magnetic Therapy:** The Meridius AcuMagnet Patch provides targeted magnetic therapy (1100 Gauss) with high-quality magnets, relieving pain and stimulating acupressure points.
- **Durable and Comfortable Design:** These water-resistant, latex-free patches are made for long-lasting use and adapt to different body areas.

Meridius AcuMagnet Patch – H type



Code	Size
MEMP03H	14 cm x 14.5 cm

Meridius AcuMagnet Patch – Coin type




Code	Size
MEMP01C	2,7 cm ø

Meridius AcuMagnet Patch – Grid type



Code	Size
MEMP04G	2.6 cm x 2 cm

Meridius AcuMagnet Patch – Rectangular type



Code	Size
MEMP02R	13 cm x 6 cm

Replacement Plasters



- Small (13 mm) is used for all kinds of press pellet
- Medium (19 mm) is used for the magnetic press needles
- Large (30mm) was created for use with the 6,000 gauss magnetic (MEA) electrodes or rubber GSR electrodes

Code	Diameter	Contents
IT5007A	13 mm	96pcs / pack
IT5007M	19 mm	120pcs / pack
IT5007B	20 mm	80pcs / pack

Magnetic Needle 3000 Gauss



- Magnetic needle, able to reach most of the auricular and body points
- The other end is rounded and may also be used for acu-pressure and acu-massage

Code	Type	Size
BN52500	3000 gauss	9 x 1.1 cm

Tsumo Shin Probe



- Spring-Activated Cutaneous Probe: Allows for precise stimulation of trigger points, acupoints, and ear points without puncturing the skin.
- Adjustable Pressure: Comes with five interchangeable springs to customize the pressure for different needs.

Code
IT6000

Spring Loaded Dermal Roller with Probe



- Multi-functional, dermal roller with spring active probe
- Rolling heads with pyramid-shaped surface

Code	Size
DB502B	Handle: 50 x 13 mm Roller: 25 x 13 mm

Regular Dermal Roller



- Sleek and simple design.
- Easy-to-hold handle with grippy texture.
- Rolling Head with tiny pyramidal-shaped surface

Code	Size
DB507-A	Handle length: 85mm
DB507-B	Handle length: 40mm

Su Jok Probes

Sharp Su Jok Probe



- Made of Chromium-nickel steel

Code	Size
DBA5302	100 x 1.5 x 2.5 mm

Round Su Jok Probe



- Made of Chromium-nickel steel

Code	Size
DBA5301	100 x 2.5 x 4.5 mm

Straight Ear Probe



Code	Size
BN5101	15cm x 2mm

Angled Ear Probe



Code	Size
BN5102	15cm x 2mm

Long Single Probe

- Very long, single probe for easy grip and application
- The handle is chiselled and offers good grip



Code

DB501A

Dual Probe

- The Dual-Probe includes pyramid spikes in the centre of the shaft.



Code

DB501B

Spring Loaded Probe

- Popular spring activated cutaneous probe
- Ideal for treating or locating trigger points



Code

DB502A

Tsun Meter Point Maker

- Interprets 1 Tsun into different distances to easily locate acupuncture points



Code

MM5001

Finger Massage Rings



Stainless steel and elastic enable the massage of fingers. Since a projection of the entire body is on each finger, the massage with a massage ring has a tangible strengthening and healing effect.

Code	Size	Colour
DBA5307BL	27 x 0.5 mm	Blue
DBA5307GR	27 x 0.5 mm	Green
DBA5307OR	27 x 0.5 mm	Orange
DBA5307RD	27 x 0.5 mm	Red
DBA5307YW	27 x 0.5 mm	Yellow
DBA5307WT	27 x 0.5 mm	White
DBA5307GD	27 x 0.5 mm	Gold

Massage Balls with Finger Ring



The massage balls are made of polystyrene. Massages hands by rolling/squeezing in a fist. Each set comes with 2 massage rings.

Code	Size	Colour
DBA5306BL	35 mm	Blue
DBA5306GR	35 mm	Green
DBA5306RD	35 mm	Red
DBA5306YL	35 mm	Yellow

Hand Massage Roller



Massage bracelet made from a resilient stainless steel wire. Used for massaging the whole hand and the wrist. The bracelets are elastic and suitable for hands of any sizes.

Code	Colour
DBA5308BL	Blue
DBA5308GR	Green
DBA5308OR	Orange
DBA5308RD	Red
DBA5308YW	Yellow
DBA5308WT	White
DBA5308GD	Gold
DBA5308SL	Silver

Small Spiral Massager



Stainless steel wire, used for the massage of both sides of the hands. The spirals are flexible and follow the hand lines, reaching even the most challenging areas.

Code	Size
DBA5305S	100 x 14 mm
DBA5305M	80 x 8.5 mm

Palm-Sized Three Channel Electric Stimulator



- For busy clinics that require several units
- Palm-sized and lightweight
- Three output channels
- Low/High intensity settings
- Unique asymmetric bi-phasic pulse
- Includes 3 sets of crocodile, needle connecting wires

Features

- High performance and versatility
- Low/High intensity settings
- Frequency selector offers a complete frequency range 1 to 500 Hz
- Suitable for needles or magnets
- Unique wave form
- Unique wave pulse minimises excessive charge build up in the tissue and damaging the nerve
- Each of the 3 channels are independently controlled
- Special dial controls allow easy use yet safe operation
- Built in battery life check with indicator light

Code

IT7002

Six-Channel Electroacupuncture Unit



- Power source: Four alkaline 'LR14' cells
- Channels: Six independently controllable
- Pulse shape: symmetric, bi-phasic square pulse
- Frequency: 0.3-500Hz
- Phase duration: 50-400 μ s, adjustable
- Amplitude: Stimulation by 'search and stimulate' probe
- Electro-acupuncture: High: 0-32mA plus/minus 25%
Low: 0-16mA plus/minus 25%
- Timer: 0-60minutes
- Output modes: (1)Constant
(2)Burst
(3)Surge
(4)Fast+Slow
(5)Sweep
(6)Random prog.1
(7)Random prog.2
(8)Random prog.3

Code

IT7001

Point Locator and Stimulator



- Four output channels
- New IC circuitry, no cross over
- Polarity reversal switch for each channel
- Sensitive point location with both light and sound indication
- Loe-Needle-Stirn treatment and location switch
- High/ Low voltage switch: low for needles, high for TENS
- Continuous, intermittent, and dense-disperse (modulated) wave forms
- Two frequency range setting: 0 to 100 and 10 to 999 Hz
- Improved alligator clip wires and detection probe
- Low battery indicator

Code
AWQ104E

PRO Acupuncture Stimulator



- Five outputs (including an output exclusively for the use of search/stimulation probe for point detection and immediate stimulation.)
- Frequency and Pulse Width are adjustable and separated into two groups, Group 1 for Ch 1 and 2, and Group 2 for Ch 3 and 4
- Two large LCD screens for displaying precise frequency of each group, and one LCD screen for showing sensitivity during acu-point detection
- Three different modes: Continuous, Modulate and Burst.
- High & Low selection of output. X1 and X5 (multiplier) of pulse rate and width adjustments
- Timer selection for different period of treatment session: Continuous, 15 minutes and 30 minutes
- Buzzer ON/OFF, for sound control during point detection
- Steady output intensity increment control, no cross over, and no interference between each output
- Music for timer alert and error warning

Code
AWQ105

Pointer Plus Locator and Stimulator



This offers sensitive acu-point or trigger point location and stimulation, which are independently adjusted with two separate control dials. The pulse frequency is fixed at 10 Hz. A flashing light and sound indicate when point detection is made.

The spring-activated detection probe allows a uniform pressure when searching. The grounding circuit is made by holding the metal of the unit or with the hand grounding pole.

- Frequency: 10Hz (fixed)
- Volts and current: 0-25 volts, 0 to 50 milliamps
- Pulse width: 220 microseconds
- Pulse shape: asymmetric biphasic square wave
- Pulse mode: continuous

Code
IT7004

Pointer Pulse



Acupuncture points locator and stimulator with infrared laser. The device measures the treatment location through sound as well as a green light indicator on the top of the unit. Location sensitivity and sound level are adjustable.

- Output channels: one
- Volts and current: Normal dial 22.5 volts, 0 to 45 milliamps Blue zone 11.0 volts, 0 to 2 microamps
- Pulse width: 260 microseconds (μ S)
- Pulse shape: asymmetric biphasic square wave with negative spike
- Pulse modes: continuous
- Laser wavelength: 650nm
- Maximum laser power: 5 mW
- Skin monitoring: sound and green light indication
- Laser operation: yellow light indication
- Power source: 9 volt alkaline battery

Code
IT7003

Pointer Plus Excel II



- Features a built-in LCD display for precise location of acupuncture and trigger points, with adjustable sensitivity and sound indicators for ease of use.
- Offers adjustable frequency (1–16 Hz) and a polarity reversal switch for tailored tonification or sedation treatments.

Code
IT7005

AS SUPER 4 digital, Needle Stimulator



- 30 fixed programs and 7 customizable programs for tailored treatments.
- 4 channels supporting up to 8 needles simultaneously.
- Adjustable intensity in 0.1 mA steps for precision and safety.
- Operates on battery or mains supply for flexibility.

Code
SM2ASS4D

Nerve and Muscle Stimulator- SDZ-II



- Independent muscle and nerve therapy.
- Suitable for various acupuncture techniques.
- Adjustable intensity and modes.
- Compact and lightweight.

Code
SMSDZ26

Pointer Pal Locator



- Easy to operate hand-held point locator. Pointer locating only, no stimulation function
- When an active point is located, a buzzing sound is produced and the pilot light on top of the unit blinks continuously

Code

IT7006

Yuwell Nerve and Muscle Stimulator



- Its small size (112cm) and portable design gives you the freedom to travel with it anywhere
- Get up to 200+ hours use when fully charged.
- Offers 10 adjustable intensity levels, catering to different needs, from mild discomfort to severe pain.

Code

SMSDP330

Pointer-Mate Acu-Locator



- The Point Mate is a pocket portable locator
- Convenient palm-size, 2.3" x 4" x 1 "
- Features sensitivity control and point location with sound and green LED light
- Includes hand ground pole and search probe with 3.5 mm jack and 9 volt battery

Code

IT7007

meridius® Manina Moxa



Meridius® Manina is a smokeless Korean moxa made from 100% real mugwort, also known as "artemesia vulgaris" or "ssuk" in Korean.

It can help reduce pain, muscle spasms and swelling and can be useful to improve blood and lymphatic circulation.

The average temperature on the skin where Meridius® Manina Moxa is applied is 43.5 degrees Celsius.

Meridius® Manina Moxa is a safe, patient-friendly product as the ash is designed to not fall off the cone. Even more, this moxa includes 40 pre-application plasters and lasts two times longer compared to standard stick-on moxa.

Code	Pieces
MA200MID	40



Stick-On Moxa

meridius® Mini Moxa



Code	Type	Contents
ME207	Smokeless	180 pcs / box



Code	Type	Contents
ME220	Natural	180 pcs / box

Meridius® Mini Moxa is a traditional, all-natural compressed type of Korean moxa shaped into small tubes. It comes in two versions: smokeless and natural. The first one produces almost no smoke or odour which makes it ideal for modern-day moxa treatments, usually performed in clinics or indoor areas.

Chosei-Kyu Stick-On Moxa

- Heat regulating paper is enclosed
- Only the highest quality ibuki moxa is used



Code	Contents
IT2001W	Weak, 100pcs / box
IT2001M	Medium, 100pcs / box
IT2001S	Strong, 100pcs / box

Neo Acu-Point Moxa

- A new type of indirect moxa
- A gentle heat penetrates acupoint



Code	Contents
IT2003A	120pcs / box
IT2003B	600pcs / box

DongBang Charcoal Moxa Cone



- The unique oven baking technique produces moxa charcoal which eliminates smoke, odour associated with traditional moxa
- The cone is coated with a special compound comprised of Bio-clay, Jade, Germanium and Bio-ceramic for enhanced effect
- It requires approximately 5 minutes until effective temperature is reached
- Highly recommended to be used with DB253

Code	Contents
DB250	approx. 55pcs / container

DongBang Moxa Device



- Screw type device for moxa treatment
- It has mesh in the middle of the device to hold the burning moxa, the heat passes up and down through the mesh-protected top and bottom
- The screw in the middle can adjust the height of the device controlling level of heat

Code
DB253

Health Preserving Moxa

- Big sized stick-on moxa.
- Moxa rolled into a cylinder with an adhesive base



Code	Contents
DB204	25 X 25 mm, 30pcs / bo

Smokeless Moxa

Chosei-Kyu Moxa 200 pcs.



- Generates no odor and little smoke
- During use, no ventilation is required, and curtains and clothes are free of odor
- A little smoke at first, but will disappear once lit
- Stable heat temperature: besides moxa, charcoal is also used as heat source

Code	Contents
IT2002S	200pcs / box


DongBang Smokeless Gu Gwan Moxa



- Excellent quality moxa for needle warming
- Smokeless and odor-free
- Charcoaled moxa rolled into a tube with hole in the middle
- It is to be used with the metal cap (DB206B Needle Cap).
- Insert the metal cap into the tube, place the metal cap on top of the acupuncture needle
- It can also be used with moxa devices (DB253 DONGBANG Moxa Device) and moxa burning bowls (DB209A and DB209B)
- Comes with 10 metal caps (DB206B)

Code	Contents
DB208A	200 pcs / box

Ondan Smokeless Needle Moxa



- Japanese, odour-free moxa
- Smokeless moxa
- 12mm long cylinder shape with a 3mm hole in the middle to go through the needle handle. Cylinder diameter is approx 12mm

Code	Contents
DBA2899-120	120pcs, 1 scoop & 6 baskets
DBA2899-480	480 cones

DongBang Gu Gwan Moxa



- Smokeless: ideal for use in a professional environment
- Odour-Free: perfect for shared clinics and sensitive patients
- High-Quality Ingredient: compressed charcoal and mugwort
- Hollow Centre: hollow for easy metal cap placement

Code	Contents
DB208B	80 pcs

Moxa Cap

- Needle cap to hold loose moxa
- Place the cap on the needle's handle



Code
DBA261

Metal Cap for Gu Gwan Moxa

- Insert the cap into the needle moxa
- Not suitable for Chinese-style needles



Code	Contents
DB206B	100 pcs / pack

Direct Super Mini Moxa

- Pre-rolled, fine thread for direct moxa treatment
- Place the burning thread directly on the skin



Code	Contents
DB212W	0.5 mm, 450pcs / box
DB212M	0.8 mm, 350pcs / box
DB212S	1.2 mm, 300pcs / box

Wakakusa Moxa – Japanese Moxa

- High-quality Japanese moxa without stems
- Great as needle moxa thanks to easy adhesion



Code	Contents
IT2021-3G	30 g
IT2021-100G	100 g
IT2021-300G	300 g

Moxa & Heating Devices

Direct & Loose Moxa

DongBang Super Gold Moxa



The best Korean pure gold moxa

Code	Contents
DB200A	10g
DB200B	50g / bag
DB200C	100g / bag

meridius® Shinhwa Gold Moxa



Meridius® Shinhwa Mogusa Gold Moxa is a loose, premium type of Korean moxa ideally used with moxa devices. This moxa is made from 100% natural mugwort.

This herb is believed to have traditional medicinal properties that could help blood circulation and increase the general state of well-being.

Code	Contents
MEDM10	10g / box
MEDM40	40g / box

Fukuju Pure Moxa Set



Medium grade moxa for everyday professional and home use. Made from good quality ibuki moxa.

Code	Contents
IT2006	10g / box

Japanese Ibuki Gold Mountain Moxa Set



- Super pure Ibuki moxa from Japan
- Easy to mould into any shape and size for direct moxa use
- Includes 5 incense sticks

Code	Contents
IT2005	10g / box

Chinese Needle Moxa



These moxa rolls can be attached securely to acupuncture needles, designed for warm-needling techniques.

Code	Contents
BN206B	200 rolls/box

Yomei Ibuki Moxa



A needle-specific moxa, Yomei moxa features super cohesiveness for easy adhesion to the needle. As a direct moxa, it is of standard quality just below that of Fukuji Wakakusa moxa.

Code	Contents
DB2031-50	50g
DB2031-100	100g

Moxa Rolls

DongBang Moxa Rolls

- Longer moxa roll with super pure moxa
- Suitable for indirect moxa treatment



Code	Contents
DB205	9 x 120 mm, 20pcs / box

Kamaya Mini Smokeless Moxa

- Japanese stick-on moxa
- The adhesive on the base helps secure the moxa



Code	Contents
DB2903-120	120pcs / box
DB2903-600	600pcs / box

Hoist Smokeless Moxa Rolls

- Very condensed charcoaled moxa
- Contains a variety of herbs



Code	Contents
BN2001	5 pcs/ box

Smokeless Moxa Perforated Rolls

- Charcoaled moxa cigar with a central hole for easy lighting.



Code	Contents
BN2020	5 pcs / box

Pure Moxa Rolls



- Classical Chinese pure moxa rolls
- For indirect moxa treatment
- Long and thick rolls with pure moxa contents

Code	Contents
DBA264	10 pcs / box

Tai-Yi Moxa Rolls



- Contains mainly moxa
- Mixed with a variety of herbs (Rhizoma Et Radiz Notopterygii, Radix Aucklandiae, Radix Angelicae Dahuricae, Radix Salviae)

Code	Contents
DBA265	10 pcs / box

Lying Smokeless Moxa Sticks



- Smokeless Moxa Roll Stick
- Suitable for most skin types
- Package included 5 pieces


Code
BN2030

Moxa & Heating Devices

Gu Gwan Moxa

meridius® Gu Gwan Moxa

- Compressed moxa formed into small tube
- Features a hollow centre in the middle



Code	Contents
ME208	7 x 25 mm, 200 pcs / box

Moxa heating device, 120 x 10 mm

- Easy to use for warming on small areas
- Warm acupoints along the meridians



Code
DB215A-J005

Lion Warmer Device

- For powerful heat application on larger areas
- Control the heat by pressing the spring on the top



Code
DB215C

meridius® Tiger Warmer Moxa

- Compatible with **Moxa Heating Device**
- Comes smokeless and natural



Code	Contents
DB215B	Natural, 5 x 80 mm, 50 pcs / box
DBA2144	Smokeless, 5 x 80 mm, 50 pcs / box

Moxa Accessories

Gu Gwan Bowl

- Ideal for indirect moxa treatment
- Features holes for acupuncture needles



Code	Type
DB209A	1 hole bowl
DB209B	5 holes bowl

OS Moxa Balm

- Can be used for direct moxa, infrared heating treatment and massage



Code	Size
MM2005	20g / jar

Shiunko Cream

- Use for direct moxa techniques, or with the tiger warmer or lion warmer



Code	Size
ITD2010S	9g / jar
ITD2010C	20g / jar

Moxa Rolls Extinguisher

- Chrome plated moxa extinguisher
- Heavy and solid



Code	Height
BN2008	40 mm

Moxa Scoop



- Ideal to remove the moxa from the needle handle
- With enclosed walls on three sides

Code	Size
BN2004	75 x 20 x 20 mm

Moxa Rolls Extinguisher with Plate



- Moxa extinguisher mounted in the middle of the plate
- Safe and stable to extinguish, to hold the moxa rolls and collect

Code	Diameter
BN2007	80 mm

Moxa Spoon



- Stainless steel spoon with a split in the middle
- Safe and easy to remove moxa from acupuncture needles

Code	Size
BN2011	Standard
BN2011S	Small

Moxa Torch



- Ideal to light moxa due to intense jet-flame
- Refillable with lighter gas (not fluid) available from most supermarkets
- Ergonomically designed for ease of holding
- Dispatched unfilled

Code
DBA257

Moxa Lighter



- Ideal to light moxa due to intense jet-flame
- Refillable with lighter gas (not fluid) available from most supermarkets
- Ergonomically designed for ease of holding
- Dispatched unfilled

Code
DBA258

Morning Star Japanese Incense



Japanese incense Nippon Koda in delicate natural fragrances

Code	Fragrance
DBA212AR	Amber, 50 sticks
DBA212CW	Cedarwood, 50 sticks
DBA212CM	Cinnamon, 50 sticks
DBA212GT	Green Tea, 50 sticks
DBA212JA	Jasmine, 50 sticks
DBA212LA	Lavander, 50 sticks
DBA212MU	Musk, 50 sticks
DBA212PA	Patchouli, 50 sticks
DBA212PN	Pine, 50 sticks
DBA212RO	Rose, 50 sticks
DBA212SW	Sandalwood, 50 sticks
DBA212VN	Vanilla, 50 sticks
DBA212YZ	Yuzu, 50 sticks
DBA212IR	Iris, 50 sticks
DBA212LO	Lotus, 50 sticks
DBA212FG	Fig, 50 sticks
DBA212GA	Gardenia, 50 sticks
DBA212MM	Mimosa, 50 sticks



Maeil Koran Incense



Sandlewood scented incence

Code	Type	Contents
DB210A	Smokeless	250 pcs / box
DB210B	Natural	200 pcs / box

Heating Lamps

TDP Infrared Heating Lamps



TDP is the abbreviation of 'the special electromagnetic spectrum' in Chinese. It has been frequently used in clinics and is also referred to as the 'Miracle Lamp'.

The TDP lamp has a plate coated with 33 important minerals for the human body. This charged thermal plate gives a spectrum between 2-50 Micro in wavelength, and between 28-34 mw/ cm2 intensity. It also gives heat element resulting in a special form of electromagnetic waves, in which both wavelength and strength are ideal for the human body.

As a result, the lamp works like moxa, giving a heating affect to depth of 3-5 cm and it is found to be very effective for patients. It can be used immediately following acupuncture and without needles, giving a similar result to the use of warm needles.

Code	Contents
HT2012UK	Table-top with digital timer
HT2000-165	Replacement plate, 16.5cm
HT2000-125	Replacement plate, 12.5cm
HT2036UK	Single Head TDP Lamp
HT2033UK	Double Head TDP Lamp

Treatment Room Equipment

Meridian & Muscle Model

- Half of the model displays meridian and points, the other half shows the muscle



Code	Size
DBA6021	80 cm
DBA6022	65 cm

Meridian Model

- Soft white vinyl with 360 points and 14 meridians
- Both Chinese and English coded points illustration



Code	Type
DBA601	Male, 50 cm
DBA602	Female, 48 cm

Hand Model

- Plastic model shows number coded acupuncture points with a booklet



Code	Size
DBA605A	22 cm
DBA605B	13 cm

Foot Model

- Made of soft vinyl
- Illustrates the acupuncture points



Code	Size
DBA607A	15 cm

Foot Model Pair of Left and Right



- A pair of soft vinyl feet that shows reflex zones of the foot
- Ideal for reflexology

Code	Size
DBA605F2	18 x 8 cm

Ear Model



- Made of soft vinyl
- Rests on a plastic base
- Both Chinese and English coded illustration
- Coding booklet include

Code	Size
DBA604A	22 cm
DBA6048	13 cm

Head Model



- Made of soft vinyl
- Features a plastic base
- Both Chinese and English coded illustration
- With a booklet
- Height 22cm

Code
DBA605

Realistic Facial Model – Silicone



- Mimics the feel of human skin, enhancing the realism of practice sessions.
- Made of high-quality silicone, ensuring long-lasting use.
- Allows practitioners to refine their skills without risk to real patients.

Code
ME153J

Animal Models



- Plastic animal model mounted on a plastic stand
- Front side of model shows acupuncture points numbered
- Comes with a booklet describing point locations and recommendations for treatments

Code	Size	Type
DBA608	26 x 7 x 23 cm	Horse
DBA609	30 x 9 x 28 cm	Dog
DBA610	28 x 8 x 19 cm	Cat
DBA612	26 x 8 x 14 cm	Cow

Dr Life Intermittent Pneumatic Compression Machines



These Intermittent Pneumatic Compression (IPC) machines are made by the world-renowned South Korean firm, Daesung Maref. Established in 1986, Daesung Maref now provide 90% of the IPC sales in the South Korean market and export to 30 countries worldwide.

With Intermittent Pneumatic Compression, the limb is surrounded by an inflatable cuff with several chambers, which are inflated sequentially so that there is a wave of compression travelling up the limb in short bursts. Venous emptying is made easy and an active circulation is expected. It is also useful in the arm to reduce lymphoedema after mastectomy.

Benefits

- Compact and simple design for home use
- Restful and calming
- Helps relieve stress for elderly and hospitalized patients
- Waist cuff can help to restore normal bowel functioning
- Helps to prevent DVT, varicose veins, lymphoedema, numbness of the extremities, poor venous return

References

- Arcelus JI et al, Home use of impulse compression of the foot and compression stocking in the treatment of Chronic venous insufficiency J Vase Surg 2007 ;34:805-11
- Kalodiki E et al, Intermittant Pneumatic Compression (IPC) in the Treatment of Peripheral Arterial Occlusive Disease (PAOD) - A Useful Tool or Just Another Device Eur J Vase Endovasc Surg 2007; 33, 309-37
- O van Bemmelen PS et al, rapid Intermittent compression increases skin circulation in chronically ischaemic legs with infra-popliteal arterial obstruction. Vasa 2000; 29: 47-52
- Labropoulos N et al, Haemodynamic effects of intermittent pneumatic compression in patients with critical limb ischaemia. J Vase Surg 2005; 42: 710 - 716

Code	Type
LX7	Dr. Life Machine
DL2003	Extention zipper

Box of 9 Standard White Couch Rolls



- 9 rolls per case
- Full 20"/508mm wide roll
- Premium quality 2ply tissue
- Easy tear perfect perforations
- Bright white in colour
- 32m per roll
- 100% pure paper
- High absorbency
- Enhanced texture
- Increased client comfort

Code
DBA6057709

Meridius Pre Injection Alcohol Swabs



- Non-woven fabric swab saturated with 70% isopropyl alcohol for skin cleansing prior to administering an injection
- Sterile and individually wrapped
- 100 pcs per box

Code	Contents
BNAS100ME	100 swabs / box

Sharps Disposal Box



- Material: Polypropylene
- Ideal for travelling, easily fitted into a bag or glove compartment
- Non-return valve

Code	Size
DBA60510-05	0.5 litre
DBA60510-10	1 litre
DBA60510-25	2.5 litre
DBA60510-375	3.75 litres
DBA60510-5	5 litres



Stainless Steel Needle tray - Small & XL



- Resistant to corrosion, rust, and stains.
- Easy to sterilize and maintain clean surfaces.
- Spill-resistant edges prevent needles from rolling off.

Code	Size
DBA5008S	110 x 35 x 10mm
DBA5008XL	140 x 65 x 15mm

Effervescent Chlorine Tablets



- Fast Acting, Fast Dissolving , Easy to Use
- Safe alternative to liquid bleach
- Lasting up to 24h, allowing you to clean as many products you wish

Code	Size
CT200	200 pcs / Box

Orders may be placed via online, e-mail, phone, or fax.

Online: www.acuprime.com

Email: info@acuprime.com

Phone: +44 (0)1392 829500

Address: 1 Forrest Units, Hennock Road East, Exeter, Devon EX2 8RU, UK,

Registered in England No. 4264020

VAT No. GB 792 9779 48

Terms & Conditions

All products shown are subject to availability.

We reserve the right to refuse to supply to any individual or company.

E-catalogue can be downloaded from www.acuprime.com.



AcuPRIME®

EXCELLENCE AT SENSIBLE PRICES™

DongBang AcuPrime Ltd

1 Forrest Units, Hennock Road East
Exeter, Ex2 8RU, United Kingdom

www.acuprime.com
info@acuprime.com

+44(0)1392 829 500